

High resolution X-ray microcalorimeter detectors with multilayer absorbers and multilayer transition-edge sensors

Mark Cunningham

Simon Labov

Matthias Frank


Troy Barbee Jr.

Lawrence Livermore National Laboratory



Microcalorimeter Approach

- Au/Bi multilayers are used to absorb X-rays
- Mo/Cu multilayer TESs are used to measure the energy of the X-rays





Mo/Cu multilayer TES

- Mo and Cu are mutually insoluble
 - Mo/Cu can be heated to 400 C to relieve stress
- Superconducting multilayers are very uniform
 - for a sharp transition, the distance between layers should be less than the length scale of the proximity effect
 - bilayers are more limited in the total thickness and resistance
 - thicker films with lower current densities are possible

Mo/Cu multilayers



Lawrence
Livermore
National
Laboratory
University of California



Total thickness = 2500 Å, 115 periods

Each period contains:

Mo: 5 monolayers @ 2.224 Å = 11.12 Å

Cu: 5 monolayers @ 2.087 Å = 10.43 Å


Thickness accuracy ± 0.2 Å
Reproducibility 0.2%

Cooper Limit Model



$$T_c = 1.14\theta_D \exp\left(\frac{d_{\text{mo}} + Bd_{\text{cu}}}{-d_{\text{mo}} * N(0)V}\right)$$

Curve fit using two variables:
 $N(0)V = 0.19$
 $B = 0.18$

Calculated $N(0)V = 0.15$



T_c versus Thickness of TES





Using multilayers, we can produce smaller pixels without increasing the current density

New Commercial SQUID Arrays




Lawrence
Livermore
National
Laboratory

University of California



36 Element 2 Stage SQUID
developed at Conductus

Noise Measurement of SQUID Array



Testing of SQUID Arrays

- Squid Arrays developed at Conductus
- No spares of sufficient quality available
from the Conductus Run
- Development continued by Robin Cantor
at STAR Cryoelectronics
- Test new arrays in a few weeks

Cu/Bi or Au/Bi Multilayer Absorber



Lawrence
Livermore
National
Laboratory
University of California

- Advantages of Cu/Bi Absorber
 - Tune Heat Capacity
 - Maintain Good Thermal Conductivity
- Bi Target Ordered
- Work have begun work on a sputtering system that will be used for the deposition of Bi

Summary

We have measured the properties of Transition Edge Devices

- T_c is constant with total thickness of the film
- T_c versus thickness ratio agrees with Cooper Limit Model
- Transitions widths are typically narrow

We are now investigating

- Stability of T_c with time and cycling
- Transition width control with annealing of samples
- Performance of new commercial SQUID system